



Resolvex™


Positive Pressure SPE

CO-DEVELOPED INSTRUMENTS & SORBENTS
FOR THE MOST CHALLENGING MATRICES
IN SAMPLE PREPARATION



What is new with Tecan ?

Positive Pressure as a universal tool for Solid Phase Extraction (SPE)



Combination of Positive Pressure Technology and Micro-particulate SPE columns results in faster sample and more consistent liquid flow compared to other commonly used technologies, e.g. as vacuum, centrifugation, gravity.

Resolvev™ – key products for ppSPE

Positive Pressure Workstations

**RESOLVEV™
M10**



Manual processor,
very robust, uniform
pressure flow!

**RESOLVEV™
A200**



Automated processor,
universal tool for sample
prep!

Smart Consumables

CEREX®



Micro Particle based SPE
for positive pressure
processing!


NBE™



Advanced SPE featuring
workflow consolidation!



Inspired by Developments in Liquid Chromatography



...and using the Power of Positive Pressure

Resolvex™ – new positive pressure SPE line

Smart Consumables

In combination with **micro particulate** sorbents technology in CEREX/NBE columns, positive pressure works as universal tool to prep **metabolites, peptides** or **fatty acids**:

- Fast liquid flow
- High pressure applications
- Micro particle sorbent
- Clean up challenging matrices
- Small sample volumes
- NBE columns integrate several workflow steps
- NBE allow micro elution

CEREX



NBE™



The SPE columns are optimized for use in positive pressure devices by its specific design and the use of micro particulate sorbents.

Classical SPE
column/plate format for
positive pressure.

The two chamber design
to prevent liquid flow
and as reaction vessel.



NBE – the all-in-one ppSPE column

Narrow Bore Extraction (NBE) Column



Reaction Reservoir enables in Column Incubation, Hydrolysis, Precipitation

Filter for initial sample clean up

Air Lock no liquid flow w/o Pressure

Microparticulate Sorbent

Small Sorbent Bed Size, high Flow Rate

Narrow bore outlet directs small elution volume

Compared to classical SPE and SLE, NBE micro particulate columns provide faster flow rates, less solvent consumption and low elution volumes.

- New two chamber design for high SPE performance
- Micro particulate sorbent for high efficiency
- Integration of pre-SPE steps in reaction vessel
- Filter for matrix clean up
- Air lock prevents dripping
- Selection of sorbents for SPE or filtration

NBE – the all-in-one ppSPE column

NBE columns – 5 steps SPE Workflow

Compared to classical SPE and SLE, micro particulate columns provide several improvements:

- faster flow rates,
- less solvent consumption
- low elution volumes.

1. Pipette sample
+ solvent

2. Protein Crash,
Incubation,
Hydrolysis

3. Apply PP to
overcome airlock
and bind analyte

4. Sorbent Bed
Drying by
Positive Pressure

5. Wash & Elute by
Positive Pressure


Positive Pressure vs. other techniques

	Positive Pressure SPE	Traditional SPE
Particles	Micro-Particulate Polymeric Sorbent ($\leq 20 \mu\text{m}$)*	Large Particle Size (30 to 70 μm)
	Mono-Disperse Particles (Typically $\pm 2 \mu\text{m}$)*	Big Particle Size Variation
	Large Contact Zone for Analyte Interaction	Solvent Channeling Through Sorbent
	Fast and Uniform	Reproducibility Effects
Pressure	Up to 5 bar Positive Pressure	Gravity, Vacuum, Centrifugation
	Fast Flow Rates	Slow Flow Rates
	Single Well Pressure Control	Vacuum Reduction
	Handles Difficult or Viscous Matrices	Clogging of Single Wells

* From Certificate of Analysis for the majority of polymeric sorbents used to manufacture T-SP SPE consumables



Available Phases



Phase Description	Phase Name
Reversed Phase	WWP, WWP2, Trace-N, Trace-N20
Mixed-Mode Anion Exchange	PSAX
Mixed-Mode Cation Exchange	PSCX
Strong Anion Exchange	HPSAX
Strong Cation Exchange	HPSCX
Weak Anion Exchange	PWAX
Weak Cation Exchange	PWCX
Special Applications	Oral Fluids: OFQC , OFXQ ; Cannabinoid: THC ; ...

Cerex and NBE consumables are optimized for use with positive pressure. The micro-particulate sorbents in plate or columns format (1, 3 & 6 ml) depending upon the phase have particle sizes typically $\leq 20 \mu\text{m}$. The sorbent materials are hydrophobic polymer-based often requiring no conditioning and equilibration steps, can be dried and are stable in organic solvents. Bed mass ranges from 2.5-50 mg depending upon column size.

Trace-N sorbent materials are silica-based C18-bonded phase (endcapped).



Application Examples

Application Area	Analyte	Matrix	Suggested Column
Vitamins	25-OH-Vitamin D	Plasma	Trace-N
Vitamins	Vitamin D Metabolites	Plasma	Maestro
Vitamins	Vitamin B	Serum Plasma	Trace-N, Maestro
Steroids	Testosterone, Cortisole and more	Serum, Plasma	WWP2
Hormones	Thyroid Hormones (rT3,T3,T4)	Serum, Plasma	THY II
Drug Monitoring	Immunosuppressant Drugs	Whole Blood	OFX
Drug Monitoring	Buprenorphine, Norbuprenorphine, and Naloxone	Urine, Oral Fluid	Trace-B, HPSCX
Drugs of Abuse, Drug Monitoring	Broad Spectrum Drugs ("Pain Panel")	Urine, Oral Fluid	OFX, HPSCX
Drugs of Abuse	Nicotine and Nicotine Metabolites	Serum, Plasma, Urine	PSCX, Trace B
Drugs of Abuse	THC, THCA, Synthetic cannabinoids	Oral Fluid, Urine	Trace-N, HPSAX, THC, PSCX
Drugs of Abuse	Amphetamines/ Metamphetamines	Urine, Whole Blood	HPSCX, Trace-B
Drugs of Abuse	Benzodiazepines	Urine, Whole Blood	HPSCX, Trace-B, Trace-J
Drugs of Abuse	EtG/EtS	Urine, Oral Fluid	EtGS, PSAX
Life Science	Peptides, metabolites or fatty acids	Serum, Plasma, Urine	NBE/CFP



Resolvex™ – new positive pressure SPE line



M10 Positive Pressure Workstation

- Positive pressure regulation from 0-5 bar
- Restricted-flow pressure manifold allows uniform pressure application
- Supports 1, 3, 6 ml (M10-48), or 1 ml cartridges, 96-well plates (M10-96)
- Operates NBE, Cerex or suitable third party columns
- Automatically adjusts to various plate heights or different columns sizes
- Column drying and slight evaporation function
- Difficult and dirty matrices
- Plug and Play, easy to operate

Fast and uniform Sample Prep with Positive Pressure

The manual M10 48/96 workstations featuring the core benefits of positive pressure technology coupled with ease of use, outstanding robustness and minimal maintenance for daily SPE sample preparation.

Resolvex™ – new positive pressure SPE line




A200 Positive Pressure Workstation

- Programmable pressure profiles from 0-5 bar
- Automated 8-channel dispenser for 11 solvents
- Restricted-flow pressure manifold allows uniform pressure application
- Supports 96 well plates or variable number of framed 1ml columns
- Operates NBE, Cerex or suitable third party columns
- Automated column drying and slight evaporation function
- For difficult and dirty matrices
- Easy to learn and freely programmable via touch screen

Fast and uniform Sample Prep with Positive Pressure

The A200 automates SPE and provides reproducible sample preparation with minimized manual intervention. The system features variable pressure profiles for gentle or strong pressure application coupled with automated dispensing of wash and elution buffers.

Key Differentiation of Tecan's Resolvex™ Technology are the microparticulate sorbents



Tecan's Resolvex™ portfolio comprises *Positive Pressure Instruments* and *Positive Pressure SPE Columns and Sorbents* that were co-developed with customers for the best possible SPE performance providing a number of clear advantages over conventional SPE.

FAQ`s for Resolvex A200

1. Since when is this device manufactured and when did it appear on the market?

- Official Launch was March 8th 2018. Manufacturing began in August 2017

2. How robust is the system? How often do breakages occur?

- The A200 is the next generation instrument following on from the previous ALD3 instrument which was on the market in the US for over 5 years. The ALD3 has proven to be a robust system and we have made design improvements to the A200 that will make it even more robust. We have also performed simulated 5 year lifetime tests on the A200.

3. Which parts of the device are most susceptible to breakdowns?

- The silicone sealing gasket on the pressure manifold may need to be replaced every 12 months depending upon usage and solvents employed. We use a Cavo syringe pump for dispensing which may need replacing every few years, depending upon usage. Besides that, we do not anticipate any other high wear components that need frequent replacing. We do recommend purchasing a service contract as you would for any automated laboratory equipment.

4. Is the interface of this device friendly to the user? Do you need some special qualification to work with it?

- The user interface is extremely simple and intuitive and uses a touch screen and or mouse. No qualification needed.

5. Which laboratories in Europe or USA already use this device? Can we get acquainted with their feedback?

- We have sold hundreds of the ALD3 but since we just launched the A200 we only have a handful of labs right now using it.

6. Are there any publications where this device is mentioned?

- We are working on clinical proteomics and steroid posters for MSACL 2018.

7. How is the process of pressure controlled, if we work with several cartridges? Can we control the pressure by installing only one SPE cartridge, with the expectation that the right pressure will then be maintained when working with multiple cartridges?

- The manifold (both 48 and 96 place) have a restrictor for each port that is controlled to allow the correct flow. The restrictor will slow flow, (increase back pressure) ensuring equal flow to all the ports regardless if there is column in that position or not. That being said it is much easier to distribute even pressure through 48 ports and the 48 place manifold have the ability to shut of unused rows with a toggle switch on the top of the manifold. The 96 place manifold will have slightly lower flow to the columns if you are only running 10 columns instead of all 96. This is because gas will flow through the path of least resistance even with the restrictor valve. This is easily compensated by increasing pressure using the regulator.



Application Support & Application Development

CORE APPLICATION TEAM EU



Dr. Christian Scherling
Business Development
Analytical Chemistry & MS



Dr. Jens Lättig
Key Account Management
Analytical Chemistry & MS



Shareef Jarvi Antar
Application Specialist
LCMS

DEMO LAB IN HAMBURG



Multiple samples (barcoded test tubes) pass through a Tecan robot → The samples are extracted into sets of multi-well plates → One set is suitable for the profiling analysis, whereas the other reformatted set is suitable for the target analysis.



Freedom EVO



Resolvex™ A200
Resolvex™ M10-96
Resolvex™ M10-48



Sciex QTRAP 5500