

Chest freezer ULT C400



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable
- Eco-friendly, HCFC-free Nature R refrigerant
- Efficient cyclopentane insulation
- ISO 9001, 14001 & 18001 OHSAS Certified

| Specifications | |
|--|-----------------------------|
| Technical | |
| Gross volume | 383 litres |
| Net volume | 368 litres |
| Temp. range | -60°C to -86°C |
| Temp. probes | 1/PT 1000 |
| Compressor | Embraco |
| Refrigerant | Nature R |
| Noise level | 55 dB |
| Power consumption 24/h | 8,10 kWh |
| Power consumption | 2.957 kWh / year |
| Voltage | 220V-240V |
| Frequency | 50 Hz |
| Features | |
| Audible & visual alarm | YES |
| Castors | 4 |
| Integrated lock | YES |
| Inner lids | 4 |
| Defrost | Manual |
| Battery Backup | 50h |
| USB | NO |
| Data connection | NO |
| Digital controller with display | XR30CX |
| Data logging | NO |
| High/Low temperature alarm | YES |
| Porthole for external probe & CO2 backup | Ø12 mm |
| Prepared for external logging equipment | YES |
| Contact for remote alarm | YES |
| Material | |
| Inner | Precoated white steel |
| Outer | Galvanized pre-coated steel |
| Dimensions | |
| Height | 850 mm |
| Height w. castors | 890 mm |
| Width | 1560 mm |
| Depth | 695 mm |
| Net weight | 88 kg |
| Internal dimensions | |
| Height | 635/385 mm |
| Width | 1190/210 mm |
| Depth | 440 mm |
| Insulation | 80 mm |
| Insulation type | PUR |
| Internal humidity | <5 % |

| ULT C400 | |
|-------------|---------|
| | |
| Product no. | N113004 |

| Storage capacity | | | | |
|------------------|--------------|--------------|--------------------|-------------|
| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
| 50 mm | 6 | 3 | 18 | N311200 |
| | 11 | 24 | 264 | N311220 |
| | | | = 282 | |
| 75 mm | 4 | 3 | 12 | N311300 |
| | 8 | 24 | 192 | N311320 |
| | | | = 204 | |
| 100 mm | 3 | 3 | 9 | N311400 |
| | 5 | 24 | 120 | N311420 |
| | | | = 129 | |