

EVO Packages 2019

VALID UNTIL DECEMBER 2019







Colony picking and Plating

EVO Packages



Colony picking and Plating EVO Packages

Platform Name	EVO75 PetriPlater 75	EVO 100 Pickolo 100	EVO 150 Pickolo 150
			
	Liquid LiHa 2, RoMa PetriPlater 75	Air LiHa 8, RoMa Pickolo Plus	
Description	A compact platform in a benchtop format to plate Petri dishes.	Colony picking and liquid handling in a compact, benchtop format.	Automate the most labor-intensive step of your cloning experiments
Application proven with	Plating of up to 60 petri plates per hour, up to 2 sample dilutions per plate.	Picks bacteria, yeast, algae, mycelia and mammalian cells...	Pick up to 576 colonies/h, template script included

For bacterial and mammalian cell colony picking

Template scripts available

Sophisticated and powerful software allows to create colony picking profiles

FLEXIBILITY

* PickoloPLUS is a SciRobotics product








Bioprocessing

EVO Packages



Bioprocessing EVO Packages

Platform	EVO 150	EVO 150	EVO 100
Name	Protein Purification RoboColumn 1	Protein Purification RoboColumn 2	Protein Purification RoboColumn 3
			
	Liquid LiHa 8, RoMa	Liquid LiHa 8, RoMa	Liquid LiHa 8
	Te-Chrom™	Te-Chrom™	Te-Chrom™
	incl. Reader!		
Description	Te-Chrom™ wizard allows easy and fast set up of protein purification processes	Te-Chrom™ wizard allows easy and fast set up of protein purification processes	Te-Chrom™ wizard allows easy and fast set up of protein purification processes
Application proven with	"Repligen"- about 200 Te-Chrom installed in the field	"Repligen"- about 200 Te-Chrom installed in the field	"Repligen"- about 200 Te-Chrom installed in the field

Create eight chromatograms in parallel

Miniaturization saves costs for resins and reagents

Take advantage of the strong partnership between Repligen and Tecan








Genomics

EVO Packages



PCR, NAP & NGS EVO Packages

Platform	EVO 100	EVO 100	EVO 100
Name	PCR	NAP	NGS
			
	Air LiHa 8 CPAC	Air LiHa 8, RoMa BioShake, Alpaqua	Air LiHa 8, RoMa ThermoShaker, Alpaqua, CPAC
Description	Effortless set-up of PCR reactions with a simple touchscreen-based interface	Wide range of tested nucleic acid purification protocols with dedicated operator interface	Verified NGS sample preparation protocols
Application proven with	Generic, not for a specific vendor protocol	Generic, not for a specific vendor protocol To learn more visit www.tecan.com/nap	To learn more visit www.tecan.com/ngs

